830 Douglas Ave. Dunedin, FL 34698 (727)733-2447 Fax:(727)733-3962 www.OceanOptics.com





Maya2000 and Maya2000Pro Data Sheet

Description

The Ocean Optics Maya2000 Series Spectrometers (Maya2000 and Maya2000Pro) includes the linear CCD-array optical bench, plus all the circuits necessary for spectrometer operation. The result is a compact, flexible system, with no moving parts, that's easily integrated as an OEM component.



The Maya spectrometers are a unique combination of technologies providing users with high sensitivity for low light-level, UV-sensitive and other scientific applications. The electronics have been designed for considerable flexibility in connecting to various modules as well as external interfaces. The Maya interfaces to PCs, PLCs and other embedded controllers through a USB 2.0 connection. The information included in this data sheet provides detailed instructions on the connection and operation of both of the Maya spectrometers.

The detector used in the Maya spectrometer is a scientific-grade, back-thinned, CCD array from Hamamatsu (product number S9840 for Maya2000 and S10420 for Maya2000Pro). For complete details on these detectors, visit <u>www.Hamamatsu.com</u>.

The Maya operates from power provided through the USB, or from a separate + 5VDC power supply. The Maya is a microcontroller-controlled spectrometer, thus all operating parameters are implemented through software interfacing to the unit.



Features

- □ Hamamatsu high UV-sensitivity detector
 - S9840 Detector for Maya2000 with Peak QE: >90%
 - S10420 Detector for Maya2000Pro with Peak QE: 75%
 - Back-thinned for good UV sensitivity
 - MPP operation for low noise operation, low dark current, wide dynamic range
- □ Spectrometer Design:
 - Symmetrical Crossed Czerny Turner
 - 101.6 mm focal length
 - 15 gratings including the HC-1 composite grating
 - 6 slit widths
- **□** Electrical Performance:
 - 16 bit, 500kHz A/D Converter
 - Integration time: 6ms 10s (Maya2000)

6ms – 5s (Maya2000Pro)

- □ Embedded microcontroller allows programmatic control of all operating parameters and standalone operation:
 - USB 2.0 480Mbps (high-speed) and 12Mbps (full speed)
 - Communication Standards for digital accessories (I2C)
- Onboard Pulse Generator:
 - 2 programmable strobe signals for triggering other devices
 - Software control of nearly all pulse parameters
- □ Onboard GPIO:
 - 10 user-programmable digital I/O
- **EEPROM** storage for:
 - Wavelength Calibration Coefficients
 - Linearity Correction Coefficients
 - Absolute Irradiance Calibration (optional)
- □ Software and Quasi Real-time triggering
- □ Plug-n-Play Interface for PC applications
- **D** 30-pin connector for interfacing to external products



Specifications

Specifications	Criteria			
Physical Specifications:				
Physical Dimensions (LxWxH)	149 mm (5.86 in.) x 109.3 mm (4.30 in.) x 50.4 mm (1.98 in.)			
Spectrometer Weight	0.96 kg (2.1 lbs.)			
Power Supply Weight	0.45 kg (1 lb.)			
Power:				
Power requirement	500 mA at +5 VDC			
Supply voltage	4.5 – 5.5 V			
Power-up time	~2s depending on code size			
Absolute Maximum Ratings:				
Vcc	+ 5.5 VDC			
Voltage on any pin	Vcc			
Spectrometer:				
Design	Symmetric crossed Czerny-Turner			
Focal length (input)	F/4 101.6 mm			
Input Fiber Connector	SMA 905 to single-strand optical fiber (0.22 NA)			
Gratings	14 different gratings			
Entrance Slit	5, 10, 25, 50, 100, or 200 μm slits. (Slits are optional. In the absence of a slit, the fiber acts as the entrance slit.)			
Detector	Maya2000: Hamamatsu S9840; Maya2000Pro: S10420			
Pixels (active)	Maya2000: 2048 x 14; Maya2000Pro: 2048 x 64			
Pixel size	14µm ²			
Spectral range	165 – 1100nm			
Quantum efficiency	Maya2000: >90% peak Maya2000Pro: 75% peak			
Well Depth	Maya2000: 130 Ke- Maya2000Pro: 200 Ke-			
Sensitivity	Maya2000: ~0.45 counts/e-; Maya2000Pro: ~0.32 counts/e-			
Dark Current	4000 e/pixel/sec (typ) @ 25° C; 200 e/pixel/sec (typ) @ 0° C			
Filters	OFLV-MAYA-200 an OFLV-MAYAPRO-200 available with HC- 1 grating			
Spectroscopic:				
Integration Time	Maya2000: 6 ms – 10s			
Dynamic Range (Typical)	Maya2000Pro: 6ms – 5s			
Dynamic Range (Guaranteed)	Maya2000: 8000:1+; Maya2000Pro: 12000:1+			
Signal-to-Noise	Maya2000: 5000:1; Maya2000Pro: 8000:1			
Dark Noise (single dark spectrum)	Maya2000: 350:1; Maya2000Pro: 450:1			
	Maya2000: 13 RMS counts; Maya2000Pro: 8.2 RMS counts			
Nonlinearity (uncorrected)	(Guaranteed)			
Linearity (corrected)	Maya2000: ~4%; Maya2000Pro: ~10%			
	>99.7%			



Maya Data Sheet

Specifications	Criteria		
Environmental Conditions: Temperature Humidity	–30° to +70° C Storage & –0° to +50° C Operation 0% – 90% noncondensing		
Interfaces: USB	USB 2.0, 480 Mbps		
(Typ. Ta=25 'C) (Typ.	(Typ. Ta=25 °C) 00 00 00 00 00 00 00 00 00 0		

Quantum Efficiency of Maya2000 S9840 Detector Quantum Efficiency of Maya2000Pro S10420 Detector

Optical Performance

The following table shows the Maya resolution for various slit sizes.

5 micron	10 micron	25 micron	50 micron	100 micron	100 micron
Slit	Slit	Slit	Slit	Slit	Slit
~1.5 pixels	~2.0 pixels	~2.5 pixels	~4.2 pixels	~8.0 pixels	~15.3 pixels

Electrical Pinout

Listed below is the pin description for the Maya Accessory Connector (J3) located on the front vertical wall of the unit. The connector is a Pak50TM model from 3M Corp. Headed Connector Part# P50-030P1-RR1-TG. Mates with part# P50-030S-TGF (requires two: 1.27mm (50 mil) flat ribbon cable: Recommended 3M 3365 Series)